


Analyte Capture from Liquid Samples: Size Matters.

Monika Weber


Comsol Conference 2013 Boston
October 10, 2013

Monika Weber, Mark A. Reed
Yale University
Department of Electrical Engineering






Current applications
of pillar structures





Toner et. al, 2007, Nature




Henderson et. al, 2006

Hwang et. al, 2011







Array of pillars




Functionalize C-R region with antibodies




Blood introduction
protein absorption



Wash with buffer
UV cleave and release



Transfer to NWs
with 2nd antibodies





Stern et. al, 2010

Lets consider a microchannel with 100 μ m pillars, injected liquid with $c_0 = 50$ nM ~ 1.6 ng/mL PSA concentration at flow rate $Q \sim 20$ nL/min very slow


The **simulated** PSA concentration after passing MPC is 45nM, which shows **10%** capture efficiency

MPC have shown **experimental 10%** capture efficiency





MPC Geometry
 Microfluidic channel
 with Array of
 functionalized
 micro-pillars



microchannel
 100um x 50um x
 10um
 Micro-pillar with
 capture sites

Mathematical Model

- Reaction at NW surface
- Convection & diff.
- Mass balance
- Navier-Stokes equation
- BC @ t=0 – instantaneous fluid injection

$$k_{ads} c + \theta \Leftrightarrow c_s$$

$$k_{des}$$

$$\frac{\partial c}{\partial t} + \nabla \cdot (-D \nabla c + \mathbf{c} \mathbf{u}) = 0$$



$$\frac{\partial c_s}{\partial t} + \nabla \cdot (-D_s \nabla c_s) = k_{ads} c (\theta_0 - c_s) - k_{des} c_s$$

$$\rho \frac{\partial \mathbf{u}}{\partial t} - \nabla \cdot \eta (\nabla \mathbf{u} + (\nabla \mathbf{u})^T) + \rho \mathbf{u} \cdot \nabla \mathbf{u} + \nabla p = 0$$


$$c = c_0 \quad \nabla \cdot \mathbf{u} = 0 \quad c_s = 0$$

Efficiency for varied structure size


Capture of same initial analyte concentration with varied pillar size



Pillar: 2 μm ; Effic. $\sim 10^{13}$




Pillar: 4 μm ; Effic. $\sim 10^5$



Pillar: 10 μm ; Effic. $\sim 10^3$

Capture on 2um pillars in the same microchannel geometry with varied initial analyte concentration



DIFFUSION dominates
HIGH capture cross section
SMALL pillar size


2um pillars

$$N_b = 2.5 \cdot 10^8$$

Diffusion 200um on that scale

Fill factor

20%




Nb - number of binding sites
Nc - number of analyte molecules

Capture efficiency varies depending on initial analyte concentration and the pillar size.

The efficiency varies from 10% to more than 99.999%.

The choice of Micro-Purification Chip geometry for sample processing has a significant influence on the accuracy of analyte concentration measurements



Size Matters

Thank you

Yale

